

# BI+BS4

## Digital Endoscope



[www.portable-ultrasound-scanner.com](http://www.portable-ultrasound-scanner.com)

**Quality, Professional, Innovative**

# Handheld Digital Medical Scopes BS4

## Specifications

### Control Unit

Resolution	1920 x 1080 pixels
Light source	Natural white (LED)
LCD screen	3.5 inch colorful TFT-LCD
Magnification	8~20 times
Image format	JPEG (photo) , H.264 (video)
Interface	Mini USB port, AV out port
Media	Media Micro SD flash card, 8G
Power supply	3.7V / 2600mAh rechargeable lithium battery
External power	AC 100~240 V , 50/60 Hz
Working time	3 hours
Charging time	5 hours
Package	305mm (L) * 260mm(W) * 65mm(H)

### Our Images




### Lens

<b>Digital Otoscope (ear inspection)</b>	
Focus range	5~30 mm
Dimensions	8.5 × 9 × 20 cm3
Weight	330g



<b>Digital Ophthalmoscope (eye inspection)</b>	
View angle	40 degree or 25 degree
Dimensions	6.5 × 9 × 20.5 cm3
Weight	340g



<b>Digital Laryngoscope (throat inspection)</b>	
Focus range	80 angle diagonal
Dimensions	8.5 × 5 × 20 cm3
Weight	310g



<b>Digital Dermatoscope (skin inspection)</b>	
Focus range	10 mm diameter
Dimensions	8.5 × 8 × 20 cm3
Weight	370g



**Wuxi Biomedical Technology Co., Ltd.**

[www.portable-ultrasound-scanner.com](http://www.portable-ultrasound-scanner.com)

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: sales@cnbimedical.com

Mobile phone: 86-13382885910

**BI+** is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.